
Instructions for Use

Volumetric infusion pump · Model PX-240

PX-240 · TECHNICAL FLOW DIAGRAM



1 Intended use

The pump delivers continuous or intermittent intravenous infusion under clinical supervision.

2 Warnings

Do not operate the device if the housing, keypad or infusion set is damaged. Use by trained clinical staff only.

3 Operating data

Keep flow, time and accuracy values in the designated table fields.

Parameter	Specification	STATUS
Flow range	0.1–1 200 mL/h	OK
Delivery accuracy	±2.5 %	OK
Pressure alarm	15–120 kPa	MON
Battery life	≥ 8 h	OK